

CONCENTRATE Wexicide 128

**A one-step cleaner disinfectant concentrate
for healthcare, schools and industry**

Germicidal Detergent
Cleans • Disinfects • Deodorizes

- A powerful one-step cleaner, disinfectant and deodorizer for decontamination of hard non-porous surfaces that have the potential to be contaminated by blood or other potentially infectious material (OPIM).
- Its excellent detergency and germ killing power especially suited to the cleaning needs of Operating Rooms, Intensive Care, Nurseries, Emergency Areas, Patient Rooms and Rest Rooms.
- Low pH makes it ideal for one-product cleaning of rest rooms and shower rooms.

HOW TO USE THE ONE PRODUCT SYSTEM

1. Add 1 ounce Wex-cide 128 to one gallon water in plastic pail.
2. Use soft sponge or cloth and the Wex-cide 128 solution to clean in the following order: Mirrors, Glass and Chrome Shelves and Dispensers; Sinks and attached Chrome Fixtures and Piping Cover; Seat and Exterior of the Toilets and Urinals; attached Fixtures and Piping; Partitions and Wall Surfaces. Surfaces must remain wet for 10 minutes.
3. Using a johnny mop and Wex-cide 128 solution, clean the interior surfaces of the toilets and urinals. To disinfect toilet bowls, add 1 pump full (1 oz. Wex-cide 128 concentrate) to the water in the toilet bowl. Allow 10 minutes before flushing.
4. Mix a solution of 1 ounce of Wex-cide 128 to each gallon of water in a mop bucket and mop all floor areas.
5. For cleaning and decontaminating surfaces soiled with blood or other body fluids which may contain HIV-1 (AIDS Virus) - see label for special instructions.

DO NOT MIX WITH OTHER CLEANERS

**GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL,
FUNGICIDAL and VIRUCIDAL in soft or hard water
(AOAC 400 ppm Hardness as Calcium Carbonate),
and organic soils (5% serum).**

WexfordLabs
The Disinfectant Solution

325 Leffingwell Avenue, Kirkwood, Missouri 63122
www.wexfordlabs.com





**A ONE-STEP CLEANER DISINFECTANT CONCENTRATE
EPA REGISTRATION NUMBER 34810-31**

GERMICIDAL QUALIFICATIONS: (AOAC Use Dilution Confirmation)

Staphylococcus aureus 1:128 (1 oz. per gallon of water)
Salmonella choleraesuis 1:128 (1 oz. per gallon of water)
Pseudomonas aeruginosa 1:128 (1 oz. per gallon of water)

ADDITIONAL ORGANISMS that passed the AOAC use dilution confirmation in 400 ppm hard water and 5% blood serum at a dilution of 1:128:

Streptococcus pyrogenes Group A and B, *Pneumococcus pneumonia*, *Shigella dysenteriae*, *Klebsiella pneumonia*, Methicillin resistant *Staphylococcus aureus* (MRSA), Methicillin resistant *Staphylococcus epidermidis* (MRSE), Tetracycline, Novobiocin, Streptomycin resistant *Staphylococcus aureus* (TNSR), Gentamicin, Methacillin resistant *Staphylococcus aureus* (GRM), Penicillin resistant *Staphylococcus* (PR), *Enterococcus Faecalis* (VRE).

TUBERCULOCIDAL ACTIVITY TEST:

AOAC, *Mycobacterium Tuberculosis* 1:128 (1 oz. per gallon of water)

FUNGICIDAL TEST:

AOAC, *Trichophyton Mentagrophytes* 1:128 (1 oz. per gallon of water)

VIRUCIDAL ACTIVITY TEST:

Influenza A3 (Hong Kong)..... 1:128 (1 oz. per gallon of water)
Adenovirus Type 2..... 1:128 (1 oz. per gallon of water)
Herpes Simplex, Types 1 and 2 1:128 (1 oz. per gallon of water)
Vaccinia Virus 1:128 (1 oz. per gallon of water)
Human Immunodeficiency Virus, Type 1 (AIDS)..... 1:128 (1 oz. per gallon of water)
Human Rotavirus 1:128 (1 oz. per gallon of water)
Respiratory Syncytial Virus 1:128 (1 oz. per gallon of water)

HARD WATER EFFICACY: All of the above tests were performed in hard water (400 ppm as Calcium Carbonate) and organic soils (5% blood serum).

EYE IRRITATION LEVEL: The concentrate is an eye irritant.

SKIN IRRITATION LEVEL: Not a primary skin irritant as defined by the Federal Hazardous Substances Act (FHSA).

ORAL TOXICITY: Not a toxic substance as defined by the Federal Hazardous Substances Act (FHSA).

DETERGENT TYPE: Synthetic-Biodegradable

APPEARANCE: Clear transparent liquid

CLEANING EFFICIENCY: Excellent

pH OF CONCENTRATE: 2.2 ± 0.4

pH OF USE-DILUTION (1:128 in deionized water): 2.8 ± 0.4